


# Actolind<sup>®</sup> w Solution

Antiseptic irrigation solution for wound  
cleansing, moisturising and care

## IT ACCELERATES THE HEALING OF WOUNDS AND BURNS



**Polyhexanide**  **Poloxamer**

A **Synergistic Approach** in Wound Care with  
**Polyhexanide + Poloxamer**

*Your best choice for wound care*

- It accelerates the wound healing process by keeping the area where it is applied clean and moist
- It is suitable for cleaning, moisturizing and decontaminating wounds and medical equipment such as bandages, dressings and wound fillers used in wound care
- Polyhexanide is classified as a first-choice active substance for the care of chronic wounds and burns colonized and infected
- It is highly effective against multi-drug resistant microorganisms (**MDRO, MRSA, VRE**) and mycotic pathogens (*Candida albicans*)
- Dermatologically tested and certified as "excellent"
- It effectively cleans scabs, necrotic tissues and biofilm layers from wounds
- It provides atraumatic change of dressing material adhering to the wound
- It prevents the re-infection of the debrided and cleaned wound
- Actolind formula is patented all over the world
- It is biocompatible and used safely
- It supports granulation and epithelization
- It can be used safely for a long time, and well tolerated by allergy and diabetes patients



Take care of the wounds with Actolind<sup>®</sup>



## Actolind® w Solution

**Actolind® w Solution** is a cleansing, irrigation, decontamination, moistening and care solution used together with materials such as bandages, dressings, pads and wound fillers for acute, chronic and infected wounds as well as 1<sup>st</sup> and 2<sup>nd</sup> degree burns.

**Actolind® w Solution** cleans wounds mechanically. Surface active agents in its composition provide a convenient and painless cleaning needed before wound care by decreasing the surface tension. **Actolind® w Solution** contributes to the wound healing process by ensuring a highly positive cleansing and moisturizing effect on applied areas.

**Actolind® w Solution** is effective against multidrug-resistant (e. g. MRSA, VRE) and mycotic pathogens. It is also well tolerated by allergics and diabetics.

### Chemical Composition:

Poloxamer, Polyhexamethylene biguanide (Polyhexanide, 0.1%) purified water

### Indications:

- Infected or uninfected acute wounds (e.g. traumatic wounds that occur as the result of cuts, tears, excoriation, crushing, being bitten, being stung or surgical wounds)
- Infected or uninfected chronic wounds (diabetic ulcers, venous ulcers, decubitus (pressure sore), etc.)
- Infected or uninfected thermal or chemical wounds (1<sup>st</sup> and 2<sup>nd</sup> degree burns, chemical burns, wounds that caused by electrical current, radiation wounds etc.)
- Entry regions of urological catheters,
- Maintenance of peristomal skin, fistulas and abscesses, wound surface prior to dressing, suitable for antiseptic cleansing, softening and removal of dressing materials, bandages, dressings, pads and wound fillers that are dried, contaminated and adhered to the wound

**Actolind® w Solution** should be used for wound surface cleaning and irrigation before any treatment. **Actolind® w Solution** can be applied in combination with **Actolind® w Gel**. The wound should be cleaned with **Actolind® w Solution** before the application of **Actolind® w Gel** to remove any contaminants and deposits.

**Actolind® w Solution** can also be used for antiseptic cleansing and softening of encrusted dressing materials that are dried and adhered to the wound surface and difficult to remove.

**Actolind® w Solution** does not inhibit granulation or epithelialization. It has neither irritating nor pain giving properties. It causes no sensitization and can be used safely for a long time.

### Usage:

**Actolind® w Solution** can be applied directly onto the wound as well as onto the bandage and gauzes.

The entire wound and its surroundings should be thoroughly cleaned with **Actolind® w Solution**. For optimal wound cleaning, it is recommended to wet the wound consisting of deep and intensively necrotic tissues, for at least 10-15 min. with the help of a fully wetted compress (application of wet dressing). The application can be repeated according to the wound conditions without any time limits. Do not dilute the solution! The cleaned wound can be covered with **Actolind® w Gel** and bandaged properly.

Thoroughly wet dressings and bandages adhered to the wound surface with **Actolind® w Solution**, wait 5-10 min. for softening and remove them carefully from wound. Wetting the dressing materials with **Actolind® w Solution** prevents adhering to the wound surface and facilitates dressing change significantly. In order to maintain a moist wound climate, the dressing can be changed daily, if necessary, several times a day. **Actolind® w Solution** is compatible with all modern wound dressings.

Possible infections can be avoided by regular cleaning catheter inlet ports, flexible connectors and its environment by **Actolind® w Solution**. Spray **Actolind® w Solution** on the area from a short distance and wait 1-5 min. The ingredients of **Actolind® w Solution** are particularly mucous membrane-friendly, cleanse catheter and PEG material without damaging it. As with wound care, encrusted sticking plasters can be removed easily after wetting with **Actolind® w Solution**.

**Actolind® w Solution** can be heated to body temperature prior to use.

Packaging				
Art. No	Unit	Packaging Type	Packaging Unit	PZN
06.6703.05	10x5 ml	Monodose Strip	100 x 5	16889754
06.6703.10	5x10 ml	Monodose Strip	50 x 5	16684727
06.6703.50	50 ml	Bottle with spray cap	100	10985385
06.6703.100	100 ml	Bottle with spray cap	81	10985391
06.6703.250	250 ml	Bottle with spray cap	49	10985416
06.6703.350	350 ml	Bottle	30	16684733
06.6703.500	500 ml	Bottle	30	10985422
06.6703.1	1000 ml	Şişe	10	10985439

### Do not use Actolind® w Solution in the following cases:

1. It is not suitable for use in hyaline cartilage and aseptic joint surgeries, in the middle and internal ear and in the eye (If **Actolind® w Solution** contacts the aseptic cartilage, wash immediately with Ringer solution or physiological saline solution),
2. If the patient is known or suspected to be hypersensitive to any of the ingredients,
3. On central nervous system or meninx,
4. Along with anionic surfactants, cleansing soaps, creams, oils, enzymes, etc. These substances should be removed fully from the wound surface before use,
5. For peritoneal lavage and washing,
6. It should not be used in the first four months of pregnancy. Depending on the risk-benefit evaluation, its use later in pregnancy should be evaluated according to the doctor's recommendation.

**This product is a Class III Medical Device according to the Medical Device Directive 93/42/EEC**

ACTO GmbH

Schulstrasse 8 • 31613 Wietzen-Holte / Germany  
Tel: +49 (0) 5022 9801-0 • Fax: +49 (0) 5022 9801-37

 [www.actogmbh.com](http://www.actogmbh.com)  
 [info@actogmbh.com](mailto:info@actogmbh.com)

