

care for your wound

Diabetic Foot Ulcer
Care and Treatment Management



Wietzen-Holte, Germany



Polyhexanide + Poloxamer

Wound cleaning, irrigation, moisturizing,
protection and care products



ACTO GmbH
Head Office & Factory
Schulstrasse 8 • 31613 Wietzen-Holte / Germany
Tel: +49 (0) 5022 9801-0 • Fax: +49 (0) 5022 9801-37
Braunschweig Office
Büchnerstrasse 11 • 38118 Braunschweig / Germany
Bayern Branch
Balthasar - Neumann Str. 4 / 97204 Höchberg / Germany
Baden Württemberg Branch
Untere Hauptstr. 6 • 78573 Wurmlingen / Germany
Tel: +49 (0) 531-239 508-0 • Fax: +49 (0) 531-239 508-11
www.actogmbh.com • info@actogmbh.com

German Company



www.actogmbh.com



CARE PRODUCTS FOR DIABETIC FOOT WOUNDS

Wound Irrigation

Wound Care Solution

- ❖ Actolind® w Solution
- ❖ Actolind® Solution Plus

✓ It provides cleaning of infected or non-infected diabetic wounds
✓ It prevents biofilm formation.
✓ It is highly effective against resistant microorganisms (MRSA, VRE, MDRO) and mycotic pathogens (*Candida albicans*)



Aclind® Sterile Gauze Compress

✓ It is used for wound cleaning of infected or non-infected diabetic ulcers.



Wound Care

Wound Care Gel

- ❖ Actolind® w Gel
- ❖ Actolind® Gel Plus



✓ It provides care for infected or non-infected diabetic ulcer wounds.
✓ It prevents biofilm formation
✓ It is highly effective against resistant microorganisms (MRSA, VRE, MDRO) and mycotic pathogens (*Candida albicans*)


Primary Wound Dressing

- ❖ Actolind® Foam Non Adhesive
- ❖ Aclind® Wound Foam
- ❖ Actolind® Plus Foam



✓ It is used as a wound dressing for partial and full thickness diabetic foot ulcers. The outermost layer is a breathable, bacteria-proof, virus-proof and waterproof polyurethane film that protects the wound from external contamination. Once the exudate is absorbed, the foam matrix forms a soft gel which allows easy removal of the bandage.

- ❖ Actolind® Hydrocolloid/Thin
- ❖ Actolind® Plus Hydrocolloid
- ❖ Actolind® Hydrocolloid Plus
- ❖ Actolind® Hydrocolloid Foam






✓ It is used as a wound dressing for partial and full thickness diabetic foot ulcers. The moist wound environment created under the dressing accelerates the natural healing process, reduces pain and protects against bacterial and other external contaminants.

Secondary Dressing

- ❖ Aclind® Sterile Gauze Compress
- ❖ Aclind® Sterile Bandage



WAGNER CLASSIFICATION AND TREATMENT

	Description	Treatment
 <div>GRADE 1</div>	Superficial ulcer	Office debridement and contact casting
 <div>GRADE 2</div>	Deeper, full thickness extension	Operative formal debridement and contact casting
 <div>GRADE 3</div>	Deep abscess formation or osteomyelitis	Operative formal debridement and contact casting

Wound Washing and Preparation of the Wound Bed

Wound Care Solution

- ❖ Actolind® w Solution
- ❖ Actolind® Solution Plus

✓ Cleaning, irrigation, moistening and care solution for diabetic foot wounds.



Wound Care Gel

- ❖ Actolind® w Gel
- ❖ Actolind® Gel Plus

✓ Protective moisturizing and care gel for diabetic foot wounds.



Aclind® Sterile Gauze Compress

✓ It is used for cleaning diabetic foot wound



Primary Dressing

Partial and full thickness diabetic foot ulcers

- ❖ Actolind® Foam Non Adhesive



- ❖ Actolind® Plus Foam



- ❖ Aclind® Wound Foam



- ❖ Actolind® Hydrocolloid Foam



- ❖ Actolind® Hydrocolloid/Thin



- ❖ Actolind® Plus Hydrocolloid



- ❖ Actolind® Hydrocolloid Plus



Secondary Dressing

- ❖ Aclind® Sterile Gauze Compress

✓ It is used in situations requiring light pressure and to all kinds of dressing materials on the wound.



- ❖ Aclind® Sterile Bandage



* Daily Dressing